

Standard Thermal Cycling: 25  $\mu$ L Endpoint PCR

- For all components except CleanAmp™ dNTPs and DNA polymerase, thaw the reaction components, vortex to mix, centrifuge briefly and store on ice.
- Prepare CleanAmp™ dNTPs:
  - Thaw at room temperature or on ice.
  - Vortex and pulse centrifuge to thoroughly mix.
  - If necessary, remove an aliquot of the stock solution and dilute with water or buffer (pH 8-10.5) to desired working concentration.
- Prepare a mastermix containing all components except for the DNA template sample. Add each of the components as shown in Table 1 (multiply amounts by the number of reactions needed) in a centrifuge tube, on ice.
- Mix the mastermix gently to protect the enzyme, by pipetting up and down. (Do not Vortex.) Pulse spin if necessary.
- Aliquot 20  $\mu$ L of mastermix into each thin-walled PCR tube.
- To each 20  $\mu$ L aliquot of mastermix, add 5  $\mu$ L of the appropriate template DNA for a final reaction volume of 25  $\mu$ L.
- Pulse spin PCR tubes to remove bubbles and collect reaction solution at bottom of tube.
- Place the tubes into a thermal cycler with a heated lid and perform the appropriate cycling conditions for target amplification:  
95°C for 10 min  
[95°C for 40 sec; X°C<sup>1</sup> for 1-30 sec; 72°C for 0.5-2 min]  
35-40 cycles 72°C for 10 min
- Analyze an aliquot of the completed reaction by agarose gel electrophoresis.

Table 1

Component	Final Concentration (25 $\mu$ L reaction)	Volume per reaction
Forward/Reverse Primer	50-500 nM	Variable
Sterile De-ionized Water	Up to 25 $\mu$ L	Up to 25 $\mu$ L
MgCl <sub>2</sub> (50 mM) <sup>2</sup>	4.5 mM	2.25 $\mu$ L
10X PCR Buffer <sup>3</sup>	1X	2.5 $\mu$ L
CleanAmp™ dUTP Mix <sup>4</sup> (10 mM)	0.2 mM	1 $\mu$ L
Taq DNA Polymerase (5 units/ $\mu$ L)	0.05 units/ $\mu$ L	0.25 $\mu$ L
DNA Template	Variable	Variable
<b>Total Volume (<math>\mu</math>L)</b>	<b>25 <math>\mu</math>L</b>	<b>25 <math>\mu</math>L</b>

<sup>1</sup> The annealing temperature should be chosen for optimal PCR performance. Most primer design software recommends an annealing temperature. The annealing temperature can also be optimized experimentally by using a thermal cycler with gradient functionality or by performing sequential experiments in which the annealing temperature is varied.

<sup>2</sup> For protocols utilizing alternate DNA polymerases, the MgCl<sub>2</sub> concentration can be optimized between 3.0 and 5.0 mM.

<sup>3</sup> 10X PCR buffer for Invitrogen's Taq DNA polymerase (200 mM Tris-HCl, pH 8.4 at 25°C; 500 mM KCl)

<sup>4</sup> CleanAmp™ dUTP should be used in a 3-fold excess over dATP, dCTP and dGTP.