

EPO mRNA encodes the human erythropoietin (EPO) protein, a hormone that controls erythropoiesis, or red blood cell production. EPO also plays a role in wound healing and the brain's response to neural injury. Karikó et al. showed that transfection of EPO mRNA *in vivo* resulted in significant increases of both reticulocyte counts and hematocrits.

This mRNA is capped using CleanCap™, TriLink's proprietary co-transcriptional capping method, which results in the naturally occurring Cap 1 structure with high capping efficiency. It is polyadenylated, modified with 5-methoxyuridine and optimized for mammalian systems. It mimics a fully processed mature mRNA.

L-7209-20 (20 µgrams)  
L-7209-100 (100 µgrams)  
L-7209-1000 (1 mg)  
L-7209-BK (Bulk amount)

1.0 mg/mL in 1mM Sodium Citrate (pH 6.4)  
mRNA Length: 858 nucleotides

Store at or below -40°C

### QC Analysis

Identity and Purity  
Agarose Gel Mobility; Pass  
Concentration: ± 6%; Pass

Product released by Quality Assurance

<sup>1</sup>A standard conversion factor of 40 µg/OD<sub>260</sub> was used to calculate quantity.

## Handling

Store at or below -40°C. Thaw and work with mRNA on ice. Upon first use, pulse spin before opening and aliquot into single use portions. Do not vortex. Use only certified RNase-free reagents and consumables with proper RNase-free technique. Use of barrier tips is recommended. Avoid freeze/thaw cycles. Do not mix with media containing serum unless first complexed with a stabilizing transfection reagent.

## References

Karikó K, Muramatsu H, Keller JM, Weissman D. Increased erythropoiesis in mice injected with submicrogram quantities of pseudouridine-containing mRNA encoding erythropoietin. *Mol Ther.* 2012 May;20(5):948-53.

Products containing the CleanCap technology are for research use only. Not for use in diagnostic or therapeutic procedures. Use of CleanCap technology may be covered by one or more patents or pending Patent Applications. [www.trilinkbiotech.com/cleancap/license.asp](http://www.trilinkbiotech.com/cleancap/license.asp).

Not for resale without express written permission. Not for use in humans. No license under any patent or patent pending is granted or implied by the purchase of any TriLink product. TriLink does not warrant that the use or sale of the products delivered hereunder will not infringe the claims of any United States or other patents or patents pending covering the use of the product alone or in combination with other products or in the operation of any process. All and any use of TriLink product is the purchaser's sole responsibility.

The Modified Nucleic Acid **Experts**™